Code No. 10777

Anti-Human Gd-IgA1(KM55) Rat IgG MoAb

- ① Paraffin-embedded sections of 3µm thickness
- 2 Deparaffinization with a series of xylene and ethanol concentrations and subsequent rehydration

Xylene 10min. x2
100% Ethanol 5min. x3
90% Ethanol 5min. x1
70% Ethanol 5min. x1
50% Ethanol 5min. x1
Rinse with distilled water 5min

- 3 Antigen retrieval: 0.05%bacterial protease subtilisin A (Sigma-Aldrich, Tokyo, Japan) dissolved in 5mmol/l tris (hydroxymethl)-aminomethane buffer (pH7.6) room temp. 2hr *This process is extremely important. If the dose of bacterial protease subtilisin A is not enough, KM55 mAb staining will not be succeeded.
- 4) Rinse with distilled water.
- ⑤ Block with Protein Block (DAKO, Tokyo, Japan) . room temp. 30min
- ⑥ 1st Ab: Anti-Human Gd-IgA1(KM55) Rat IgG 100μg/mL 37°C 60min
- Several washes with phosphate-buffered saline and triethanolamine-buffered saline containing 0.05% tween-20.
- ② 2nd Ab: Alexa Fluor 555(or Other colors)-conjugated goat anti-rat IgG antibody (1:1000; Life Technologies, Carlsbad, CA). 37°C 30min
- Several washes with phosphate-buffered saline and triethanolamine-buffered saline containing 0.05% tween-20.
- Seal with Fluoromount (Diagnostic BioSystems, Pleasanton, CA).