

ELISA

The ELISA is one of the most important and popular method to detect and measure cytokines. The assay is fast, easy to perform and is one of the most sensitive and reliable technologies available.

U-CyTech offers high-quality ELISA kits and antibody pairs for the detection of cytokine and granzyme levels in biological fluids, such as cell culture supernatant, plasma or serum.

The kits contain coating antibodies, biotinylated detection antibodies, cytokine standards, streptavidin-HRP conjugate, TMB substrate solution, stop solution, cytokine stabilization buffer, BSA, Tween-20 and ELISA plates with cover slips.

